INSTRUCTIONS FOR THE USE OF YOUR NIGHT GUARD

- It may take a little time for you to get used to your nightguard, but we recommend that you wear it every night.
- After each use, rinse your nightguard with cool water. You may use a toothbrush and liquid soap to clean it several times a week, but do not use toothpaste to clean it. It is too abrasive.
 - -When your nightguard is not in use, it should be stored dry in its case.

If you have any questions please feel free to call our office.

John E. Smith, D.M.D. Michael J. Welch, D.D.S. 38 S Last Chance Gulch Helena, MT 59601 406-443-6160